



Implant Details: Height Definition

Catalog #	Implant Size (Length x Post. Height x Lordosis)	Posterior Height (mm)
5203-20424L6	24 mm x 4 mm x 6°	4
5203-20624L6	24 mm x 6 mm x 6°	6
5203-20824L6	24 mm x 8 mm x 6°	8
5203-21024L6	24 mm x 10 mm x 6°	10
5203-21224L6	24 mm x 12 mm x 6°	12
5203-20428L6	28 mm x 4 mm x 6°	4
5203-20628L6	28 mm x 6 mm x 6°	6
5203-20828L6	28 mm x 8 mm x 6°	8
5203-21028L6	28 mm x 10 mm x 6°	10
5203-21228L6	28 mm x 12 mm x 6°	12
5203-20424L12	24 mm x 4 mm x 12°	4
5203-20624L12	24 mm x 6 mm x 12°	6
5203-20824L12	24 mm x 8 mm x 12°	8
5203-21024L12	24 mm x 10 mm x 12°	10
5203-20428L12	28 mm x 4 mm x 12°	4
5203-20628L12	28 mm x 6 mm x 12°	6
5203-20828L12	28 mm x 8 mm x 12°	8
5203-21028L12	28 mm x 10 mm x 12°	10
5203-20424L18	24 mm x 4 mm x 18°	4
5203-20624L18	24 mm x 6 mm x 18°	6
5203-20824L18	24 mm x 8 mm x 18°	8
5203-20428L18	28 mm x 4 mm x 18°	4
5203-20628L18	28 mm x 6 mm x 18°	6
5203-20828L18	28 mm x 8 mm x 18°	8

Note: Heights are measured from Tooth Root to Tooth Ro 0.5 mm tall, adding approximately 1 mm to the overall heig